








ABRAXIS® 2,4-D Coated Tube ELISA 54004B

<p>1.</p>  <p>Add 500 µL of either Standards (calibrators), Control or Samples to the bottom of each test tube.</p>	<p>5.</p>  <p>Add 500 µL of Color Reagent down the inside wall of each test tube.</p>
<p>2.</p>  <p>Add 200 µL of 2,4-D Enzyme Conjugate to each test tube.</p> <p>Vortex or swirl.</p>	<p>6.</p>  <p>Incubate for 20 minutes.</p> <p>Prepare blank.</p>
<p>3.</p>  <p>Incubate for 20 minutes.</p>	<p>7.</p>  <p>Add 500 µL of Stopping Solution to each test tube.</p> <p>Read OD 450</p>
<p>4.</p>  <p>Flood tubes with deionized water.</p> <p>Invert tubes and blot.</p> <p>Repeat this step 4 times.</p>	