

Ivermectin in Honey Sample Preparation

1. Intended Use

For the detection of Ivermectin in honey

2. Materials and Reagents Required

Deionized or distilled water

Acetonitrile

50 mL Centrifuge tube

Centrifuge

Methanol

Micropipettes with disposable plastic tips/Serological pipettes

Vortex mixer

ABRAXIS® Ivermectin Plate ELISA Kit (PN 5142B)

3. Sample Preparation Procedure

3.1 Heat up honey sample to make it homogenous.

3.2 Weigh 5 g of honey in 50 mL centrifuge tube.

3.3 Add 10 mL of 80% Acetonitrile/Water in the tube.

3.4 Vortex for 2 minutes.

3.5 Shake for 3 minutes by hand for complete extraction of Ivermectin.

3.6 Centrifuge for 5 minutes at 6000 rpm.

3.7 Transfer 1 mL of the supernatant and dry under gentle stream of nitrogen gas until complete dry.

3.8 Add 1 mL of 60% MeOH/Water and reconstitute the dried extract.

3.9 Employ 50 µL per well in ABRAXIS® Ivermectin Plate ELISA Kit.

* Use dilution factor of **1.6** to calculate the concentration of Ivermectin.

4. For ordering or technical assistance contact:

Gold Standard Diagnostics

795 Horsham Road

Horsham, PA 19044

WEB: www.abraxiskits.com

Phone: (215) 357 3911

Fax: (215) 357 5232

Ordering: info.abraxiskits@us.goldstandarddiagnostics.com

Technical Support: support.abraxiskits@us.goldstandarddiagnostics.com

Date this Technical Bulletin is effective: 05/16/2024

Version: 01