

Saxitoxin in Lobster Tomalley Sample Preparation

1. Intended Use

For the detection of Saxitoxin in lobster tomalley.

2. Materials and Reagents Required

Strainer (#10 sieve)

Plastic tablecloth (for protecting work area)

Deionized or distilled water (for rinsing samples prior to homogenization)

Immersion blender or equivalent electronic blender

600 mL plastic beaker

Permanent marker

10% Bleach solution (for cleaning equipment between samples)

Centrifuge

ABRAXIS[®] Saxitoxin (PSP), ELISA, 96-test (PN 52255B)

3. Notes and Precautions

This procedure is intended for use with lobster tomalley. Other matrices should be thoroughly validated before use with this procedure.

4. Sample Preparation and Extraction Procedure

NOTE: If a 100 g sample is needed for regulatory purposes, extraction solution volume should be adjusted accordingly

4.1 Remove tomalley from the lobster, wash with deionized water and homogenize.

4.2 Mix 10 g of homogenized tomalley with 10 mL of 0.1M HCl and boil for 5 minutes while stirring.

4.3 Allow to cool. Centrifuge for 10 minutes at approximately 3500 g.

4.4 Collect supernatant. Adjust pH to < pH 4.0 with 5 N HCl.

4.5 Remove 10 µL and dilute in 10 mL of 1X Sample Diluent (this will be a 1:1,000 dilution). Vortex

4.6 Analyze as sample (Assay Procedure, step 1 in the user guide)

5. Evaluation of Results

The ABRAXIS[®] Saxitoxin (PSP) ELISA Kit can be evaluated using commercial ELISA evaluation programs, and the concentrations of samples are determined using a standard curve of various Saxitoxin concentrations run with each test. Detailed information on the evaluation of results using the ABRAXIS[®] Saxitoxin (PSP) ELISA Kit can be found in the ABRAXIS[®] Saxitoxin (PSP) ELISA Kit User Guide.

6. For ordering or technical assistance contact

Gold Standard Diagnostics

795 Horsham Road

Horsham, PA 19044

WEB: www.abraxiskits.com

Phone: (215) 357 3911

Fax: (215) 357 5232

Ordering: info.abraxaxis@us.goldstandarddiagnostics.com

Technical Support: support.abraxaxis@us.goldstandarddiagnostics.com

Date this Technical Bulletin is effective: 05/16/2024

Version: 01